

ADVANCED SUTURE INSTRUMENTS

Precision in Every Stitch



Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



CRILE-MURRAY

		TC
150 mm – 6 „	07-101-150	07-102-150
180 mm – 7 "	07-101-180	
200 mm – 8 "	07-101-200	



BAUMGARTNER

		TC
120 mm – 4 ¾ „	07-103-120	07-104-120
140 mm – 5 ¾ "	07-103-140	07-104-140



Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



CRILE-WOOD



		TC	TC
150 mm – 6 "	07-105-150	07-106-150	07-106-151
180 mm – 7 "	07-105-180	07-106-180	07-106-181
200 mm – 8 "	07-105-200	07-106-205	
230 mm – 9 "		07-106-230	
260 mm – 10 ¼ "		07-106-265	
300 mm – 12 "		07-106-300	

glatt
smooth
lisse
lisa
liscio




1:2
BABY-CRILE-WOOD

	TC
150 mm – 6 "	07-107-150

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



BROWN

130 mm – 5 "	07-108-130
180 mm – 7 "	07-108-180



COLLIER

125 mm – 5 "	07-109-120
--------------	------------

kräftiges Modell, Maul diagonal gerieft
heavy pattern, diagonally serrated jaw
modèle fort, mors strié diagonalement
modelo fuerte, boca estriada diagonalmente
modello forte, morso con strie diagonali



1:2
DERF



1:2



1:2
DERF (WRIGHT)

		TC
125 mm – 5 "	120 mm – 4 ¾ "	120 mm – 4 ¾ "
07-110-120	07-111-120	07-112-125

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



WEBSTER

		TC	TC
130 mm – 5 "	07-113-130	07-114-130	07-115-130

glatt
smooth
lisse
lisa
liscio

glatt
smooth
lisse
lisa
liscio

gerieft
serrated
strié
estriada
zigrinato

NEIVERT


125 mm – 5 "	07-116-130
--------------	------------

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi




140 mm – 5 ½ "	160 mm – 6 ¼ "	180 mm – 7 "	200 mm – 8 "	240 mm – 9 ½ "	260 mm – 10 ¼ "	300 mm – 12 "
07-122-140	07-122-160	07-122-180	07-122-200	07-122-240	07-122-260	07-122-300



MAYO-HEGAR

TC	TC	TC	TC	TC	TC	TC	TC
140 mm – 5 ½ "	160 mm – 6 ¼ "	180 mm – 7 "	200 mm – 8 "	240 mm – 9 ½ "	260 mm – 10 ¼ "	300 mm – 12 "	340 mm – 13 ½ "
07-123-140	07-123-160	07-123-180	07-123-200	07-123-240	07-123-260	07-123-300	07-123-340



MAYO-HEGAR

140 mm – 5 ½ "	160 mm – 6 ¼ "	180 mm – 7 "	200 mm – 8 "	240 mm – 9 ½ "	260 mm – 10 ¼ "
07-124-140	07-124-160	07-124-180	07-124-200	07-124-240	07-124-260

schlankes Modell
delicate pattern
modèle fin
modelo delicado
modello delicato

MAYO-HEGAR

TC	TC	TC	TC	TC	TC
140 mm – 5 ½ "	160 mm – 6 ¼ "	180 mm – 7 "	200 mm – 8 "	240 mm – 9 ½ "	260 mm – 10 ¼ "
07-125-140	07-125-160	07-125-180	07-125-200	07-125-240	07-125-260

schlankes Modell
delicate pattern
modèle fin
modelo delicado
modello delicato

MAYO-HEGAR



TC	TC
240 mm – 9 ½ "	260 mm – 10 ¼ "
07-126-240	07-126-260

gebogen
curved
courbe
curva
curvo



Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



1:2
HEGAR

1:2

		TC
150 mm – 6 "	07-127-150	
175 mm – 7 "	07-127-175	
205 mm – 8 "	07-127-205	07-128-200
240 mm – 9 1/2 "		07-128-240

kräftiges Modell
heavy pattern
modèle fort
modelo fuerte
modello forte

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



HEANEY

		TC
160 mm – 6 ¼ "	07-129-160	07-130-160
180 mm – 7 "	07-129-180	07-130-180
200 mm – 8 "	07-129-210	07-130-210

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



SAROT

		TC
180 mm – 7 "	324-316-180	07-132-180
200 mm – 8 "		07-132-200
260 mm – 10 ¼ "	324-316-260	07-132-260



JAMESON 1:2

		TC
235 mm – 9 ¼ "	07-133-230	07-134-230

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



MASSON

		TC
270 mm – 10 ¾"	07-135-270	07-136-270



1:2




FINOCHIETTO

		TC
200 mm – 8"	07-137-200	07-138-200
270 mm – 10 ¾"	07-139-270	07-140-270



Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi




1:1



1:1

STRATTE

		TC
210 mm – 8 1/4"		07-142-210
230 mm – 9"	07-141-230	07-142-230



1:2

WANGENSTEEN

1:2

		TC
270 mm – 10 3/4"	07-143-270	07-144-270

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



BOZEMANN (WERTHEIM)

		TC
200 mm – 8 "	07-145-200	07-146-200
240 mm – 9 1/2 "	07-145-240	07-146-240
260 mm – 10 1/4 "		07-146-260
300 mm – 12 "		07-146-300




1:2
THOMSON (BOZEMANN)

200 mm – 8 "	07-147-200
--------------	------------

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



GILLIES

145 mm – 5 ¾" 07-148-145

GILLIES

		TC
130 mm – 5"		07-149-130
160 mm – 6 ¼"	07-148-160	07-149-160



FOSTER-GILLIES

	TC	TC	TC	TC
120 mm – 4 ¾"	07-150-120	07-152-120	07-153-120	07-155-120



1:2



1:2



1:2



1:2

gerade, gerieft
straight, serrated
droit, strié
recta, estriada
retto, zigrinato

gerade, glatt
straight, smooth
droit, lisse
recta, lisa
retto, liscio

gebogen, gerieft
curved, serrated
courbe, strié
curva, estriada
curvo, zigrinato

gebogen, glatt
curved, smooth
courbe, lisse
curva, lisa
curvo, liscio

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



OLSEN-HEGAR

		TC
115 mm – 4 ½ "	07-156-115	07-157-120
140 mm – 5 ½ "	07-156-140	07-157-140
160 mm – 6 ¼ "	07-156-160	07-157-160
180 mm – 7 "	07-156-180	07-157-180

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



MATHIEU

140 mm – 5 ½"	07-158-140	07-159-140	07-160-140
170 mm – 6 ¾"	07-158-170	07-159-170	07-160-170
200 mm – 8"	07-158-200	07-159-200	07-160-200
240 mm – 9 ½"	07-158-240		

mit Längkehle im Maul
grooved jaw
boca acanalada
mors strié
morso scannellato

schlankes Modell
delicate pattern
modèle fin
modelo delicado
modello delicato



Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



MATHIEU-OLSEN

		TC
140 mm – 5 ½"	07-163-140	07-164-140
170 mm – 6 ¾"	07-163-170	07-164-170
200 mm – 8"	07-163-200	07-164-200



Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



CASTROVIEJO-
JACOBSON



140 mm – 5 ½ "	07-168-140	07-169-140	07-170-140	07-171-140
180 mm – 7 "	07-168-180	07-169-180	07-170-180	07-171-180
210 mm – 8 ¼ "	07-168-210	07-169-210	07-170-210	07-171-210

glattes Maul
smooth jaws
mors lisse
boca lisa
morso liscio

gerade, ohne Sperre
straight, without ratchet
droit, sans crémaillère
recta, sin cremallera
retto, senza cremagliera

gerade, mit Sperre
straight, with ratchet
droit, avec crémaillère
recta, con cremallera
retto, con cremagliera

gebogen, ohne Sperre
curved, without ratchet
courbe, sans crémaillère
curva, sin cremallera
curvo, senza cremagliera

gebogen, mit Sperre
curved, with ratchet
courbe, avec crémaillère
curva, con cremallera
curvo, con cremagliera



CASTROVIEJO-
JACOBSON



	TC	TC	TC	TC
140 mm – 5 ½ "	07-172-140	07-173-140	07-174-140	07-175-140
180 mm – 7 "	07-172-180	07-173-180	07-174-180	07-175-180
210 mm – 8 ¼ "	07-172-210	07-173-210	07-174-210	07-175-210

glattes Maul
smooth jaws
mors lisse
boca lisa
morso liscio

gerade, ohne Sperre
straight, without ratchet
droit, sans crémaillère
recta, sin cremallera
retto, senza cremagliera


gerade, mit Sperre
straight, with ratchet
droit, avec crémaillère
recta, con cremallera
retto, con cremagliera

gebogen, ohne Sperre
curved, without ratchet
courbe, sans crémaillère
curva, sin cremallera
curvo, senza cremagliera

gebogen, mit Sperre
curved, with ratchet
courbe, avec crémaillère
curva, con cremallera
curvo, con cremagliera

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



TOENNIS

		TC
180 mm – 7 "	07-176-180	07-177-180






LICHTENBERG

		TC
140 mm – 5 ½ "		07-179-140
170 mm – 6 ¾ "	07-178-170	07-179-170
200 mm – 8 "	07-178-200	07-179-200

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



Nadelhalter
Needle holders


Porte-aiguilles
Porta-agujas
Porta aghi



RYDER

	TC
130 mm – 5 "	07-192-130
150 mm – 6 "	07-192-150
180 mm – 7 "	07-192-180
200 mm – 8 "	07-192-200
230 mm – 9 "	07-192-230
260 mm – 10 ¼ "	07-192-260

MINI-RYDER



	TC
07-198-130	
07-198-150	
07-198-180	
07-198-200	
07-198-230	
07-198-260	

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



KILNER

		TC
135 mm – 5 ¼"	07-199-135	07-200-135
160 mm – 6 ¼"	07-199-160	



1:2

CONVERSE

		TC
130 mm – 5"	07-201-125	



HÖSEL

TC	TC
200 mm – 8"	250 mm – 9 ¾"
07-202-200	07-202-250

mit langem Hals
with long neck
avec col longue
con cuello largo
con collo lungo



TC	TC
205 mm – 8"	260 mm – 10 ¼"
07-203-200	07-203-260

mit langem Hals
with long neck
avec col longue
con cuello largo
con collo lungo

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi

TC



ZWEIFEL

TC	TC
200 mm – 8 "	230 mm – 9 "
07-204-200	07-204-230




LANGENBECK

	TC
140 mm – 5 ½ "	07-205-140
160 mm – 6 ¼ "	07-205-160
180 mm – 7 "	07-205-180
200 mm – 8 "	07-205-200
230 mm – 9 "	07-205-230

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



DE BAKEY

TC	TC	TC	TC	TC	TC
160 mm – 6 ¼"	180 mm – 7"	200 mm – 8"	230 mm – 9"	260 mm – 10 ¼"	300 mm – 12"
07-206-160	07-206-180	07-206-200	07-206-230	07-206-260	07-206-300




	TC	TC
150 mm – 6"	07-207-150	07-208-150
180 mm – 7"	07-207-180	07-208-180
200 mm – 8"	07-207-200	07-208-200

Drahtzange

Wire forceps
Pince pour fils d'acier
Pinza para hilos
Pinze per fili d'acciaio

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi



VASCULAR

TC	TC	TC	TC
140 mm – 5 ½ "	180 mm – 7 "	200 mm – 8 "	230 mm – 9 "
07-209-140	07-209-180	07-209-200	07-209-230



v



	TC	TC	TC	TC
180 mm – 7 "	07-210-180	07-211-180	07-212-180	07-213-180

gerade, ohne Sperre
straight, without ratchet
droit, sans crémaillère
recta, sin cremallera
retto, senza cremagliera

gerade, mit Sperre
straight, with ratchet
droit, avec crémaillère
recta, con cremallera
retto, con cremagliera

gebogen, ohne Sperre
curved, without ratchet
courbe, sans crémaillère
curva, sin cremallera
curvo, senza cremagliera

gebogen, mit Sperre
curved, with ratchet
courbe, avec crémaillère
curva, con cremallera
curvo, con cremagliera



DIETRICH

	TC
160 mm – 6 ¼ "	07-214-160
180 mm – 7 "	07-214-180
200 mm – 8 "	07-214-200
230 mm – 9 "	07-214-230
250 mm – 9 ¾ "	07-214-260

Nadelhalter
Needle holders

Porte-aiguilles
Porta-agujas
Porta aghi




CASTROVIEJO



130 mm – 5 "	19-922-130	19-923-130
--------------	------------	------------

gerade, gerieftes Maul
straight, serrated jaws
droit, mors strié
recta, boca estriada
retto, morso zigrinato



gebogen, gerieftes Maul
curved, serrated jaws
courbe, mors strié
curva, boca estriada
curvo, morso zigrinato

CASTROVIEJO

	TC	TC
130 mm – 5 "	20-724-130	20-725-130

gerade, gerieftes Maul
straight, serrated jaws
droit, mors strié
recta, boca estriada
retto, morso zigrinato

gebogen, gerieftes Maul
curved, serrated jaws
courbe, mors strié
curva, boca estriada
curvo, morso zigrinato



CASTROVIEJO


	TC
130 mm – 5 "	07-215-130

glattes Maul
smooth jaws
mors strié
boca lisa
morso liscio



BARRAQUER

	TC	TC
130 mm – 5 "	07-216-130	07-217-130



gebogen, mit Sperre
curved, with ratchet
courbe, avec crémaillère
curva, con cremallera
curvo, con cremagliera



gebogen, ohne Sperre
curved, without ratchet
courbe, sans crémaillère
curva, sin cremallera
curvo, senza cremagliera

Nahtinstrumente
Suture instruments

Instruments à suture
Instrumentos para sutura
Strumenti per sutura



REVERDIN


	Fig. 1	Fig. 2	Fig. 3	Fig. 4
150 mm – 6 "	07-232-150	07-233-150	07-234-150	07-235-150
195 mm – 7 ¾ "	07-236-001	07-236-002	07-236-003	07-236-004
225 mm – 8 ¾ "	07-232-230	07-233-230	07-234-230	07-235-230



YOUNG

250 mm – 9 ¾ "	07-237-000
----------------	------------

Komplett mit einem Nadelaufsatz
Complete with one needle
Complet avec une aiguille
Completo con una aguja
Completo con un ago



07-237-010	07-237-020	07-237-030
------------	------------	------------

Ersatznadel
Spare needle
Aiguille de rechange
Aguja de cambio
Ago di ricambio

Nahtinstrumente
Suture instruments

Instruments à suture
Instrumentos para sutura
Strumenti per sutura



200 mm – 8"

07-238-200

biegsam
flexible
malléable
maleable
flessibile



COOPER

190 mm – 7 ½"

07-239-190



KRONECKER

130 mm – 5"

07-240-130

07-241-130

für die rechte Hand
for right hand
pour la main droite
para la mano derecha
per la mano destra

für die linke Hand
for left hand
pour la main gauche
para la mano izquierda
per la mano sinistra

unterbindungsnadeln
Ligature needles

Aiguilles à suture
Agujas para sutura
Aghi per sutura



DESCHAMPS

200 mm – 8 "	07-242-200	07-243-200	07-244-200	07-245-200
240 mm – 9 ½ "	07-246-240	07-247-240	07-248-240	07-249-240
270 mm – 10 ¾ "	07-246-270	07-247-270	07-248-270	07-249-270

scharf
sharp
pointue
aguda
acuta

scharf
sharp
pointue
aguda
acuta

stumpf
blunt
mousse
roma
smussa

stumpf
blunt
mousse
roma
smussa

für die rechte Hand
for right hand
pour la main droite
para la mano derecha
per la mano destra

für die linke Hand
for left hand
pour la main gauche
para la mano izquierda
per la mano sinistra

für die rechte Hand
for right hand
pour la main droite
para la mano derecha
per la mano destra

für die linke Hand
for left hand
pour la main gauche
para la mano izquierda
per la mano sinistra

Nahtinstrumente
Suture instruments

Instruments à suture
Instrumentos para sutura
Strumenti per sutura



ROTTER



260 mm – 10 ¼"	07-250-010	07-250-020	07-251-010	07-251-020
----------------	------------	------------	------------	------------

für die rechte Hand
for right hand
pour la main droite
para la mano derecha
per la mano destra

für die linke Hand
for left hand
pour la main gauche
para la mano izquierda
per la mano sinistra

für die rechte Hand
for right hand
pour la main droite
para la mano derecha
per la mano destra


für die linke Hand
for left hand
pour la main gauche
para la mano izquierda
per la mano sinistra



SCHMIEDEN

1:2

210 mm – 8 ¼" 07-252-210



CLAUS-
EICKEN

1:2

230 mm – 9" 07-253-230




CLAUS

1:2

230 mm – 9" 07-254-230

Nahtinstrumente
Suture instruments


Instruments à suture
Instrumentos para sutura
Strumenti per sutura



HURD			DUPUY-WEISS		
1:2		1:2	1:2		1:2
210 mm – 8 ¼"	07-255-210	07-256-210	220 mm – 8 ¾"	07-257-220	07-258-220
für die rechte Hand for right hand pour la main droite para la mano derecha per la mano destra	für die linke Hand for left hand pour la main gauche para la mano izquierda per la mano sinistra		für die rechte Hand for right hand pour la main droite para la mano derecha per la mano destra	für die linke Hand for left hand pour la main gauche para la mano izquierda per la mano sinistra	

Nahtinstrumente
Suture instruments

Instruments à suture
Instrumentos para sutura
Strumenti per sutura



y



170 mm – 6 ¾"

07-259-010

07-259-020

07-259-030


07-259-040

07-259-050

07-259-060

Septumnadeln

Septum suture needles
Aiguilles de sutures du septum
Agujas de sutura de septum
Aghi di sutura per setto



07-259-070

07-259-080

07-259-090

NAGER

1:2

240 mm – 9 ½"

07-260-240



FALK


1:2

240 mm – 9 ½"

07-261-240

1:2

07-262-240



für die rechte Hand
for right hand
pour la main droite
para la mano derecha
per la mano destra


für die linke Hand
for left hand
pour la main gauche
para la mano izquierda
per la mano sinistra

führungshohlsonden – Nahtinstrumente
Ligature conductors – Suture instruments


Sondes conductrices – Instruments à suture
Sondas conductoras – Instrumental para sutura
Sonde conduttrici – Strumenti per sutura




1:3
280 mm – 11 " 4 mm 07-263-280



1:3
KIRSCHNER
255 m – 10 " 5 mm 07-264-250



1:3
BRUNNER
300 mm – 12 " 07-265-300 07-266-300



1:3
BRUNNER
8 mm 8 mm
300 mm – 12 " 07-267-300 07-268-300

führungshohlsonden – Nahtinstrumente
Ligature conductors – Suture instruments

Sondes conductrices – Instruments à suture
Sondas conductoras – Instrumental para sutura
Sonde conduttrici – Strumenti per sutura



BRUNNER


180 mm – 7 "	07-273-180
270 mm – 10 ¾ "	07-273-270

Ferrozell

1:2
KIRSCHNER


270 mm – 10 ¾ "	07-272-270
-----------------	------------

Ferrozell



KÖNIG


	2.5 mm	5 mm	7 mm
190 mm – 7 ½ "	07-269-030	07-270-190	
200 mm – 8 "			07-271-200



SCHMIEDEN

190 mm – 7 ½ "
8 mm
07-274-190

1:2




SCHMIEDEN-
PAYR

190 mm – 7 ½ "
8 mm
07-275-190

1:2

führungshohlsonden – Nahtinstrumente
Ligature conductors – Suture instruments

Sondes conductrices – Instruments à suture
Sondas conductoras – Instrumental para sutura
Sonde conduttrici – Strumenti per sutura wieN



PAYR

	5 mm	5 mm
220 mm – 8 ¾"	07-276-220	07-277-220



1:2
MOD

235 mm – 9 ¼"
07-278-230

Fadenfänger
Ligature catcher
Porte-ligature
Captador de hilos
Portalegature




1:1



1:2
BOZEMANN

210 mm – 8 ¼"
07-279-210

Schlingendreher
Knot tier
Tords-fil
Tuerce-lazos
Torcifilo



1:1



1:2

230 mm – 9"
07-280-230

Fadengabel
Ligature adjuster
Pousse-fil
Horquilla de hilos
Forchetta da sutura

1:1

klammernaht-instrumente
Suture clip instruments

Instruments pour sutures aux agrafes
Instrumentos para suturas con grapas
Strumenti per suture con agraiffe



CHILDE

1:2

180 mm – 7 "

07-281-175

kräftig, mit Klammerschiene
heavy pattern, with clip holder
modèle fort, avec porte-agrafe
modelo fuerte, con barra para agrafes
modello forte, con portaagrafe

MICHEL

1:2

120 mm – 4 ¾ "

07-282-120



1:1

1:1

WACHENFELD

1:2

125 mm – 5 "

07-283-125

07-284-125



1:1

1:2

HEGENBARTH

125 mm – 5 "

07-285-125



1:2

COLLIN

130 mm – 5 "

07-286-130



1:2

RICHTER-HEATH

140 mm – 5 ½ "

07-287-140

klammernaht-instrumente
Suture clip instruments

Instruments pour sutures aux agrafes
Instrumentos para suturas con grapas
Strumenti per suture con agraiffe



MICHEL

1:2

120 mm – 4 ¾" 07-288-120



1:1

MICHEL



1:1



1:1



1:1



7,5 mm x 1,75 mm	11 mm x 2 mm	12 mm x 3 mm	14 mm x 3 mm
07-289-007	07-289-011	07-289-012	07-289-014

Verkaufseinheit = 1.000 Stück (10 Packungen à 100 Stück)
Unit for sale = 1.000 pieces (10 packages with 100 pieces)
Unité de vente = 1.000 pièces (10 paquets avec 100 pièces)
Unidad de venta = 1.000 piezas (10 paquetes con 100 piezas)
Unità di vendita = 1.000 pezzi (10 pacchetti con 100 pezzi)



1:1

MICHEL



1:1



1:1



1:1

16 mm x 3 mm	18 mm x 3 mm	20 mm x 3 mm	22 mm x 3 mm
07-289-016	07-289-018	07-289-020	07-289-022