



# UROLOGY SURGICAL INSTRUMENTS

*Precision in Urological Care*




Nierensteinzangen  
Kidney stone forceps

Pinces pour les calculs rénaux  
Pinzas para cálculos renales  
Pinze per calcoli renali



RANDALL

| Fig. 1      | Fig. 2        | Fig. 3        |
|-------------|---------------|---------------|
| 230 mm – 9" | 220 mm – 8 ¾" | 190 mm – 7 ½" |
| 17-101-001  | 17-101-002    | 17-101-003    |



| Fig. 4        | Fig. 5        | Fig. 6        |
|---------------|---------------|---------------|
| 190 mm – 7 ½" | 190 mm – 7 ½" | 190 mm – 7 ½" |
| 17-101-004    | 17-101-005    | 17-101-006    |


Gallenstein-Instrumente

Gall stone instruments

Instruments pour calculs biliaires

Instrumental para coger cálculos biliares

Strumenti per canali bilari



**17-102-000**

Satz 9 St. 52-240-001 – 52-240-009

Set, 9 pieces 52-240-001 – 52-240-009

Jeu, 9 pièces 52-240-001 – 52-240-009

Juego, 9 piezas 52-240-001 – 52-240-009

Serie, 9 pezzi 52-240-001 – 52-240-009

**17-103-000**

Satz 13 St. 52-240-001 – 52-240-013

Set 13 pieces 52-240-001 – 52-240-013

Jeu 13 pièces 52-240-001 – 52-240-013

Juego 13 piezas 52-240-001 – 52-240-013

Serie 13 pezzi 52-240-001 – 52-240-013

Gallenstein-Instrumente

Gall stone instruments

Instruments pour calculs biliaires

Instrumental para coger cálculos biliares

Strumenti per canali bilari



DESJARDINS



|               | Fig. 1     | Fig. 2     | Fig. 3     |
|---------------|------------|------------|------------|
| 280 mm – 11 " | 17-104-001 | 17-104-002 | 17-104-003 |

federnd  
elastic  
élastique  
elástica  
elastico



|               | Fig. 1     | Fig. 2     | Fig. 3     |
|---------------|------------|------------|------------|
| 280 mm – 11 " | 17-105-001 | 17-105-002 | 17-105-003 |

federnd  
elastic  
élastique  
elástica  
elastico

Gallenstein-Instrumente

Gall stone instruments

Instruments pour calculs biliaires

Instrumental para coger cálculos biliares

Strumenti per canali bilari




DESJARDINS



|               | Fig.1      | Fig. 2     | Fig. 3     | Fig. 4     |
|---------------|------------|------------|------------|------------|
| 300 mm – 12 " | 17-106-001 | 17-106-002 | 17-106-003 | 17-106-004 |

federnd  
elastic  
élastique  
elástica  
elastico



x



1:1

MOYNIHAN

|                 |            |
|-----------------|------------|
|                 | x = 5 mm   |
| 340 mm – 13 ½ " | 17-107-340 |


biegsam  
malleable  
malléable  
maleable  
flessibile

1:3

Gallensteinlöffel  
Gall stone scoops


Curettes pour calculs biliaires  
Cucharillas para coger cálculos biliares  
Cucchiai per colcoli bilari

LUER-KOERTE



|                 | Fig. 000   | Fig. 00    | Fig. 0     | Fig. 1     |
|-----------------|------------|------------|------------|------------|
| 320 mm – 12 ½ " | 17-108-010 | 17-108-020 | 17-108-030 | 17-108-040 |

biegsam, stumpf  
malleable, blunt  
malléable, mousse  
maleable, roma  
flessibile, smussa



|                 | Fig. 2     | Fig. 3     | Fig. 4     | Fig. 5     |
|-----------------|------------|------------|------------|------------|
| 320 mm – 12 ½ " | 17-108-050 | 17-108-060 | 17-108-070 | 17-108-080 |

biegsam, stumpf  
malleable, blunt  
malléable, mousse  
maleable, roma  
flessibile, smussa



MAYO



|               |            |
|---------------|------------|
| 280 mm – 11 " | 17-109-270 |
|---------------|------------|

1:2



FERGUSON

|                |                |                |
|----------------|----------------|----------------|
| 215 mm – 8 ½ " | 220 mm – 8 ¾ " | 235 mm – 9 ¼ " |
| 17-110-215     | 17-110-220     | 17-110-235     |



1:2



Gallenstein-Instrumente

Gall stone instruments

Instruments pour calculs biliaires

Instrumental para coger cálculos biliares

Strumenti per canali bilari



BLAKE

1:1

1:1

200 mm – 8"

17-111-200

17-112-200



DESJARDINS

1:1

1:1

230 mm – 9"

17-113-001

17-113-002

Gallenstein-Instrumente

Gall stone instruments

Instruments pour calculs biliaires

Instrumental para coger cálculos biliares

Strumenti per canali bilari



CZERNY


1:1



245 mm – 9 ¾"

17-114-250


1:2



TURNER-  
WARWICK

1:1

1:1



180 mm – 7"

190 mm – 7 ½"

17-115-180

17-115-190

Prostatahaken  
**Prostatic retractors**

Ecarteurs prostatiques  
Separadores para la próstata  
Divaricatori prostatici



YOUNG

210 mm – 8 ¼" 17-116-210




YOUNG

220 mm – 8 ¾" 17-117-220



YOUNG

210 mm – 8 ¼" 17-118-210



YOUNG


220 mm – 8 ¾" 17-119-220

Instrumente für Prostata-operationen  
**Prostatectomy instruments**

Instruments pour prostatectomie  
Instrumentos para prostatectomía  
Strumenti per operazioni della prostata




MILLIN



240 mm – 9 ½" 17-120-240

**Fadenführungszange**  
Ligatur carrying forceps  
Pince porte-fil  
Pinza porta-ligaduras  
Pinza per fili di sutura



MILLIN

230 mm – 9" 17-121-230

**Fasszange für die vorderen Lappen**  
Prostatic lobe holding forceps  
Pince pour saisir les lambeaux  
Pinza para coger los lóbulos anteriores  
Pinza da presa per lobo prostatico

1:2

Instrumente für Prostata-operationen  
**Prostatectomy instruments**



Instruments pour prostatectomie  
Instrumentos para prostatectomía  
Strumenti per operazioni della prostata



**YOUNG**

200 mm – 8" 17-123-200

**Prostata-Fasszange**  
Prostate holding forceps  
Pince pour saisir la prostate  
Pinza para coger la próstata  
Pinza da presa per lobo prostatico



**MILLIN**

|             | ATRAUMA       |
|-------------|---------------|
| 230 mm – 9" | 235 mm – 9 ¼" |
| 17-124-230  | 17-125-230    |

**Kapselfasszange zur retropubischen Prostatektomie**  
Capsule forceps for retropubic prostatectomy  
Pince pour saisir les capsules pour la prostatectomie rétropubique  
Pinza capsular para la prostatectomía retropúbica  
Pinza a capsula per prostatectomia retropubica

Blasenwundspreizer  
Bladder retractors

Ecarteurs vésicaux  
Separadores vesicales  
Divaricatori vescicali



#### THOMSON-WALKER

17-126-000

komplett mit 2 Paar Seitenvalven, 1 Mittelvalve  
complete with 2 pairs of lateral blades, 1 central blade  
complet, avec 2 paires de valves latérales, 1 valve centrale  
completo con 2 pares de valvas laterales, 1 valva central  
completo, con 2 paia di valve laterali, 1 valva centrale

#### MILLIN



|                                                         |                                                                                                | 57 mm x 25 mm                                            | 80 mm x 25 mm                                            | 120 mm x 45 mm                                                                             |
|---------------------------------------------------------|------------------------------------------------------------------------------------------------|----------------------------------------------------------|----------------------------------------------------------|--------------------------------------------------------------------------------------------|
| 17-127-000                                              | 17-127-010                                                                                     | 17-127-012                                               | 17-127-014                                               | 17-127-016                                                                                 |
| komplett<br>complete<br>complet<br>completa<br>completa | Sperrer, allein<br>Spreader, only<br>Retracteur, seul<br>Retractor, sólo<br>Divaricatore, solo | 2 Blätter<br>2 Blades<br>2 Valves<br>2 Valvas<br>2 Valve | 2 Blätter<br>2 Blades<br>2 Valves<br>2 Valvas<br>2 Valve | 1 Mittelblatt<br>1 Center blade<br>1 Valve centrale<br>1 Valva central<br>1 Valva centrale |

17-128-500

komplett mit Kaltlichtaufsatz  
complete with cold light carrier  
complet avec guide de lumière froide  
completa con guía de luz fría  
completa con portaluze fredda

17-128-510

Kaltlichtaufsatz  
Cold light carrier  
Guide de lumière froide  
Guía de luz fría  
Porta luce fredda, solo

Blasenwundspreizer  
**Bladder retractors**

Ecarteurs vésicaux  
Separadores vesicales  
Divaricatori vescicali



LEGUEU

17-129-000



JUDD-MASSON

17-130-000



LEGUEU

1:3

|                  |            |
|------------------|------------|
|                  | x = 35 mm  |
| 270 mm – 10 3/4" | 17-131-270 |



KOCHER

1:3

|                  |            |
|------------------|------------|
|                  | x = 67 mm  |
| 260 mm – 10 1/4" | 17-132-270 |



MILLIN 1:3

|                  |
|------------------|
| 285 mm – 11 1/4" |
| 17-133-280       |

dilatationsbougies  
Dilating bougies  
Bougies dilatatrices  
Bujías dilatadoras  
Candelette dilatatrici

**DITTEL**

|           |                 | 220 mm – 8 ¾ " |
|-----------|-----------------|----------------|
| 2 ⅔ mm Ø  | 8 Charr. / Fr.  | 17-134-008     |
| 3 mm Ø    | 9 Charr. / Fr.  | 17-135-090     |
| 3 ⅓ mm Ø  | 10 Charr. / Fr. | 17-134-010     |
| 3 ½ mm Ø  | 11 Charr. / Fr. | 17-135-110     |
| 4 mm Ø    | 12 Charr. / Fr. | 17-134-012     |
| 4 ⅓ mm Ø  | 13 Charr./Fr.   | 17-135-130     |
| 4 ½ mm Ø  | 14 Charr. / Fr. | 17-134-014     |
| 5 mm Ø    | 15 Charr. / Fr. | 17-135-150     |
| 5 ⅓ mm Ø  | 16 Charr. / Fr. | 17-134-016     |
| 5 ½ mm Ø  | 17 Charr. / Fr. | 17-135-170     |
| 6 mm Ø    | 18 Charr. / Fr. | 17-134-018     |
| 6 ⅓ mm Ø  | 19 Charr. / Fr. | 17-135-190     |
| 6 ½ mm Ø  | 20 Charr. / Fr. | 17-134-020     |
| 7 mm Ø    | 21 Charr. / Fr. | 17-135-210     |
| 7 ⅓ mm Ø  | 22 Charr. / Fr. | 17-134-022     |
| 7 ½ mm Ø  | 23 Charr. / Fr. | 17-135-230     |
| 8 mm Ø    | 24 Charr. / Fr. | 17-134-024     |
| 8 ⅓ mm Ø  | 25 Charr. / Fr. | 17-135-250     |
| 8 ½ mm Ø  | 26 Charr. / Fr. | 17-134-026     |
| 9 mm Ø    | 27 Charr. / Fr. | 17-135-270     |
| 9 ⅓ mm Ø  | 28 Charr. / Fr. | 17-134-028     |
| 9 ½ mm Ø  | 29 Charr. / Fr. | 17-135-290     |
| 10 mm Ø   | 30 Charr. / Fr. | 17-134-030     |
| 10 ⅔ mm Ø | 32 Charr. / Fr. | 17-135-320     |



1:2

gerade  
straight  
droit  
recta  
rette

17-135-000

Satz 23 St, 8 – 30 Charr. / Fr.  
Set 23 pieces, 8 – 30 Charr. / Fr.  
Jeu 23 pièces, 8 – 30 Charr. / Fr.  
Juego 23 piezas, 8 – 30 Charr. / Fr.  
Serie 23 pezzi, 8 – 30 Charr. / Fr.

**DITTEL**

|          |                 | 220 mm – 8 ¾ " |
|----------|-----------------|----------------|
| 2 ⅔ mm Ø | 8 Charr. / Fr.  | 17-136-008     |
| 3 mm Ø   | 9 Charr. / Fr.  | 17-136-009     |
| 3 ⅓ mm Ø | 10 Charr. / Fr. | 17-136-010     |
| 3 ½ mm Ø | 11 Charr. / Fr. | 17-136-011     |
| 4 mm Ø   | 12 Charr. / Fr. | 17-136-012     |
| 4 ⅓ mm Ø | 13 Charr. / Fr. | 17-136-013     |
| 4 ½ mm Ø | 14 Charr. / Fr. | 17-136-014     |
| 5 mm Ø   | 15 Charr. / Fr. | 17-136-015     |
| 5 ⅓ mm Ø | 16 Charr. / Fr. | 17-136-016     |
| 5 ½ mm Ø | 17 Charr. / Fr. | 17-136-017     |
| 6 mm Ø   | 18 Charr. / Fr. | 17-136-018     |
| 6 ⅓ mm Ø | 19 Charr. / Fr. | 17-136-019     |
| 6 ½ mm Ø | 20 Charr. / Fr. | 17-136-020     |
| 7 mm Ø   | 21 Charr. / Fr. | 17-136-021     |
| 7 ⅓ mm Ø | 22 Charr. / Fr. | 17-136-022     |
| 7 ½ mm Ø | 23 Charr. / Fr. | 17-136-023     |
| 8 mm Ø   | 24 Charr./Fr.   | 17-136-024     |
| 8 ⅓ mm Ø | 25 Charr. / Fr. | 17-136-025     |
| 8 ½ mm Ø | 26 Charr. / Fr. | 17-136-026     |
| 9 mm Ø   | 27 Charr. / Fr. | 17-136-027     |
| 9 ⅓ mm Ø | 28 Charr. / Fr. | 17-136-028     |
| 9 ½ mm Ø | 29 Charr. / Fr. | 17-136-029     |
| 10 mm Ø  | 30 Charr. / Fr. | 17-136-030     |



1:2

gebogen  
curved  
courbe  
curva  
curve

17-136-000

Satz 23 St, 8 – 30 Charr./Fr.  
Set 23 pieces, 8 – 30 Charr. / Fr.  
Jeu 23 pièces, 8 – 30 Charr. / Fr.  
Juego 23 piezas, 8 – 30 Charr. / Fr.  
Serie 23 pezzi, 8 – 30 Charr. / Fr.



dilatationsbougies  
Dilating bougies  
Bougies dilatatrices  
Bujías dilatadoras  
Candelette dilatatrice

**GUYON**

|             |                 |            |
|-------------|-----------------|------------|
| 2 2/3 mm Ø  | 8 Charr. / Fr.  | 17-137-080 |
| 3 1/3 mm Ø  | 10 Charr. / Fr. | 17-137-100 |
| 3 3/4 mm Ø  | 11 Charr. / Fr. | 17-137-110 |
| 4 mm Ø      | 12 Charr. / Fr. | 17-137-120 |
| 4 1/3 mm Ø  | 13 Charr. / Fr. | 17-137-130 |
| 4 3/4 mm Ø  | 14 Charr. / Fr. | 17-137-140 |
| 5 mm Ø      | 15 Charr. / Fr. | 17-137-150 |
| 5 1/3 mm Ø  | 16 Charr. / Fr. | 17-137-160 |
| 5 3/4 mm Ø  | 17 Charr. / Fr. | 17-137-170 |
| 6 mm Ø      | 18 Charr. / Fr. | 17-137-180 |
| 6 1/3 mm Ø  | 19 Charr. / Fr. | 17-137-190 |
| 6 3/4 mm Ø  | 20 Charr. / Fr. | 17-137-200 |
| 7 mm Ø      | 21 Charr. / Fr. | 17-137-210 |
| 7 1/3 mm Ø  | 22 Charr. / Fr. | 17-137-220 |
| 7 3/4 mm Ø  | 23 Charr. / Fr. | 17-137-230 |
| 8 mm Ø      | 24 Charr. / Fr. | 17-137-240 |
| 8 1/3 mm Ø  | 25 Charr. / Fr. | 17-137-250 |
| 8 3/4 mm Ø  | 26 Charr. / Fr. | 17-137-260 |
| 9 mm Ø      | 27 Charr. / Fr. | 17-137-270 |
| 9 1/3 mm Ø  | 28 Charr. / Fr. | 17-137-280 |
| 9 3/4 mm Ø  | 29 Charr. / Fr. | 17-137-290 |
| 10 mm Ø     | 30 Charr. / Fr. | 17-137-300 |
| 10 1/3 mm Ø | 31 Charr. / Fr. | 17-137-310 |
| 10 3/4 mm Ø | 32 Charr. / Fr. | 17-137-320 |
| 11 mm Ø     | 33 Charr. / Fr. | 17-137-330 |
| 11 1/3 mm Ø | 34 Charr. / Fr. | 17-137-340 |
| 11 3/4 mm Ø | 35 Charr. / Fr. | 17-137-350 |
| 12 mm Ø     | 36 Charr. / Fr. | 17-137-360 |
| 12 1/3 mm Ø | 37 Charr. / Fr. | 17-137-370 |
| 12 3/4 mm Ø | 38 Charr. / Fr. | 17-137-380 |
| 13 mm Ø     | 39 Charr. / Fr. | 17-137-390 |
| 13 1/3 mm Ø | 40 Charr. / Fr. | 17-137-400 |

**GUYON**

17-134-000

Satz 27 St, 10 – 36 Charr. / Fr. ohne Innengewinde  
Set 27 pieces, 10 – 36 Charr. / Fr. without inner thread  
Jeu 27 pièces, 10 – 36 Charr. / Fr. sans filet intérieur  
Juego 27 piezas, 10 – 36 Charr. / Fr. sin filete interior  
Serie 27 pezzi, 10 – 36 Charr. / Fr. senza filetto interno

**VAN BUREN**

17-135-000

Satz 15 St, 8 – 36 Charr. / Fr.  
Set 15 pieces, 8 – 36 Charr. / Fr.  
Jeu 15 pièces, 8 – 36 Charr. / Fr.  
Juego 15 piezas, 8 – 36 Charr. / Fr.  
Serie 15 pezzi, 8 – 36 Charr. / Fr.




**VAN BUREN**

|             |                 |            |
|-------------|-----------------|------------|
| 2 2/3 mm Ø  | 8 Charr. / Fr.  | 17-135-008 |
| 3 1/3 mm Ø  | 10 Charr. / Fr. | 17-135-010 |
| 4 mm Ø      | 12 Charr. / Fr. | 17-135-012 |
| 4 1/3 mm Ø  | 14 Charr. / Fr. | 17-135-014 |
| 5 1/3 mm Ø  | 16 Charr. / Fr. | 17-135-016 |
| 6 mm Ø      | 18 Charr. / Fr. | 17-135-018 |
| 6 3/4 mm Ø  | 20 Charr. / Fr. | 17-135-020 |
| 7 1/3 mm Ø  | 22 Charr. / Fr. | 17-135-022 |
| 8 mm Ø      | 24 Charr. / Fr. | 17-135-024 |
| 8 3/4 mm Ø  | 26 Charr. / Fr. | 17-135-026 |
| 9 1/3 mm Ø  | 28 Charr. / Fr. | 17-135-028 |
| 10 mm Ø     | 30 Charr. / Fr. | 17-135-030 |
| 10 2/3 mm Ø | 32 Charr. / Fr. | 17-135-032 |
| 11 1/3 mm Ø | 34 Charr. / Fr. | 17-135-034 |
| 12 mm Ø     | 36 Charr. / Fr. | 17-135-036 |

dilatationsbougies – Metallkatheter  
Dilating bougies – Metal catheters

Bougies dilatatrices – Sondes métalliques  
Bujías dilatadoras – Sondas metálicas  
Candelette dilatatrici – Sonde metalliche



GUYON

420 mm x 175 mm x 40 mm

366-196-000

komplett, Charr. 10 – 30, im Edelstahl-Etui  
complete, Charr. 10 – 30, in steel case  
complet, Charr. 10 – 30, en étui métallique  
completa, Charr. 10 – 30, en caja metálica  
completa, Charr. 10 – 30, in astuccio metallico



|            |                 |            |
|------------|-----------------|------------|
| 2 2/3 mm Ø | 8 Charr. / Fr.  | 17-136-008 |
| 3 1/3 mm Ø | 10 Charr. / Fr. | 17-136-010 |
| 4 mm Ø     | 12 Charr. / Fr. | 17-136-012 |
| 4 2/3 mm Ø | 14 Charr. / Fr. | 17-136-014 |
| 5 1/3 mm Ø | 16 Charr. / Fr. | 17-136-016 |
| 6 mm Ø     | 18 Charr. / Fr. | 17-136-018 |
| 6 2/3 mm Ø | 20 Charr. / Fr. | 17-136-020 |
| 7 1/3 mm Ø | 22 Charr. / Fr. | 17-136-022 |
| 8 mm Ø     | 24 Charr. / Fr. | 17-136-024 |


**Metallkatheter für Frauen**  
Metal catheters for women  
Sondes métalliques pour femmes  
Sondas metálicas para mujeres  
Sonde metalliche per donne

|            |                 |            |
|------------|-----------------|------------|
| 2 2/3 mm Ø | 8 Charr. / Fr.  | 17-137-008 |
| 3 1/3 mm Ø | 10 Charr. / Fr. | 17-137-010 |
| 4 mm Ø     | 12 Charr. / Fr. | 17-137-012 |
| 4 2/3 mm Ø | 14 Charr. / Fr. | 17-137-014 |
| 5 1/3 mm Ø | 16 Charr. / Fr. | 17-137-016 |
| 6 mm Ø     | 18 Charr. / Fr. | 17-137-018 |
| 6 2/3 mm Ø | 20 Charr. / Fr. | 17-137-020 |
| 7 1/3 mm Ø | 22 Charr. / Fr. | 17-137-022 |
| 8 mm Ø     | 24 Charr. / Fr. | 17-137-024 |

**Metallkatheter für Männer**  
Metal catheters for men  
Sondes métalliques pour hommes  
Sondas metálicas para hombres  
Sonde metalliche per uomini

Urologie-Instrumente  
Instruments for urology

Instruments pour l'urologie  
Instrumentos para la urología  
Strumenti per urologia



MATHIEU

|                  |            |
|------------------|------------|
| x = 200 mm – 8"  | 17-139-200 |
| x = 280 mm – 11" | 17-139-280 |

**Fremdkörperzange**

Foreign body forceps  
Pince à corps étrangers  
Pinza para cuerpos extraños  
Pinza per corpi estranei



GUYON:3

1:3

|                |                |
|----------------|----------------|
| 390 mm – 15 ½" | 350 mm – 13 ¾" |
| 17-136-390     | 17-137-350     |

**Katheterspanner**

Mandrin for rubber catheters  
Mandrin pour cathéters  
Mandril para catéteres  
Introduttore per catheteri



1:2

|               |            |
|---------------|------------|
| 160 mm – 6 ¼" | 17-135-001 |
|---------------|------------|

**Samenstrang-Haken**

Hook for the spermatic cord  
Crochet pour le cordon spermatique  
Granchio para el cordón espermático  
Uncino per dotto spermatico




1:2

|               |            |
|---------------|------------|
| 145 mm – 5 ¾" | 17-138-125 |
|---------------|------------|

**zum Fixieren des Samenstranges**

Fixation clamp for the spermatic cord  
Pince à vis pour fixation du cordon spermatique  
Pinza con rosca para fijación el cordón espermático  
Pinza con vite per fissare il funicolo spermatico



Beschneidungsinstrumente  
Circumcision instruments

Instruments de circoncision  
Instrumentos para circuncisión  
Strumenti per circuncisione



|         |            |
|---------|------------|
| 8 mm Ø  | 17-139-008 |
| 11 mm Ø | 17-139-011 |
| 13 mm Ø | 17-139-013 |
| 16 mm Ø | 17-139-016 |
| 21 mm Ø | 17-139-021 |
| 26 mm Ø | 17-139-026 |
| 32 mm Ø | 17-139-032 |

**Beschneidungsklemme**  
Circumcision instrument  
Instrument de circoncision  
Instrumento para circuncisión  
Strumento di circuncisione



MORGAN

1:2

|            |            |            |
|------------|------------|------------|
| 75 mm – 3" | 17-140-001 | 17-140-002 |
|------------|------------|------------|

**Beschneidungsklemme**  
Circumcision instrument  
Instrument de circoncision  
Instrumento para circuncisión  
Strumento di circuncisione

**für Erwachsene**  
for adults  
pour adultes  
para adultos  
per adulti

**für Kinder**  
for children  
pour enfants  
para niños  
per bambini



17-141-103

**Zange zum Erfassen der Glans und der Vorhaut**  
Forceps for grasping the glans and the foreskin  
Pince pour saisir le prépuce  
Pinza para agarrar el prepucio  
Pinza da presa d'il prepuzio



STOCKMANN:2

17-142-000

**Penisklemme**  
Penis clamp  
Pince pour le pénis  
Clamp para el pene  
Serrapene



STRAUSS:2

17-143-000

**Penisklemme**  
Penis clamp  
Pince pour le pénis  
Clamp para el pene  
Serrapene